





ILV

**Institut Lavoisier
de Versailles**




XRD (X-RAY DIFFRACTION)



X-ray diffraction provides information on the solid state compounds (powder, crystals, thin films...). This technique is used for determining the atomic and molecular structure of crystalline materials.

Expertise



The main activity of the X-ray diffraction platform is the structural characterization of small molecules (organic molecules, coordination complexes, polyoxometalates, metal-organic frameworks)

Powder:

- » Qualitative analysis of powder : phase identification
- » Quantitative analysis of mixtures
- » Microstructurale study (size and strain parameters)
- » Pattern matching
- » Elucidation of crystal structures
- » Rietveld refinement

Single-crystals :


- » Selection and mounting
- » Determination of the unit cell parameters
- » Data collection & data treatment
- » Elucidation of the crystal structure
- » Refinement
- » Deposit of the CIF file in database

The X-ray diffraction platform is open to all research teams of our institute. Most diffraction devices are self-service after training by the manager. The platform also conducts studies for other external laboratories, public or private, and for companies. For any analysis request, please contact one of the manager.

Facilities



Single-crystal diffractometer Bruker AXS D8 Venture (2015)




- four-circles diffractometer (Kappa geometry)
- Mo microfocus source
- Photon 100 detector (CMOS technology)
- Cryostream 800+ (Liquid nitrogen ; 80K - 500K)



Single-crystal diffractometer Bruker AXS X8 APEX II (2004)

- Four-circles goniometer (Kappa geometry)
- Mo source
- APEX II 4K detector (CCD technology)
- KRYO-FLEX (Liquid nitrogen; 90K -



300K)



**Powder X-ray
diffractometer Siemens D5000
(1996)**

- Two-circles goniometer (Thêta/2
Thêta geometry)
- Cu source
- Graphite monochromator
- Scintillation detector



**Powder diffractometer Bruker AXS
D8 Advance (2004)**




- Two-circles goniometer (Thêta/2 Thêta goniometer)
- Cu source
- Ge monochromator
- LynxEye detector



Powder diffractometer Bruker AXS D8 Advance (2013)

- Thêta/Thêta geometry
- Cu source
- LynxEye detector
- HTK 1200N (25°C – 1200°C)



THE X-RAY DIFFRACTION PLATFORM IS OPEN TO THE SCIENTIFIC COMMUNITY

Manager :
Jérôme Marrot